












# eIFU Single Use Instruments



Dokument: GEMMA\_eIFU\_V02, Stand 05-2020

 <b>Manufacturer</b>	GEMMA MEDICAL Ltd. Mattenstrasse 11, CH-2555 Bruegg, Switzerland Tel.: +41 323747700 e-mail: info@gemma-medical.ch www.gemma-medical.ch
<b>Logistic</b>	Logistic, Customer Service: GEMMA MEDICAL Service GmbH, Grmuender Strasse 6, 75181 Pforzheim, Germany Tel.: +49 7231 7801215 e-mail: gemma@diatec-pforzheim.de
<b>Distributor</b>	1stQ GmbH, Konrad-Zuse-Ring 23, 68163 Mannheim, Germany Tel: +49 621 78953790 e-mail: <a href="mailto:sales@1stq.de">sales@1stq.de</a> www.1stq.de
<b>For use by</b>	Only for use by suitably trained and qualified medical professionals
<b>Intended Use</b>	The instruments are intended for ophthalmological-surgical use. Hooks, spatulas: hand-held instruments for the manipulation of ocular tissue or the artificial intraocular lens to be implanted during cataract surgery Choppers: hand-held instruments for manipulation or sculpting of ocular tissue like the natural lens Forceps: hand-held instruments to grasp or hold ocular tissue or other medical devices such as sutures or IOLs Speculums: surgical instrument for retraction of the eyelids
<b>Contraindication</b>	The single use instruments are not intended for use in surgery other than ophthalmology. The instruments are delivered sterile and may not be reprocessed or reused, see MDD statement. Use before expiration date.
<b>MDD Reuse Statement</b>	If this product is reprocessed and/or reused, GEMMA MEDICAL cannot guarantee the functionality, material structure, or cleanliness or sterility of the product. Reuse could lead to illness, infection and/or injury to the patient or user. This product is labeled as 'single use' which is defined as a device intended to be used once only for a single patient.
<b>Warning</b>	Used surgical instruments are contaminated. Disease may occur from injury. Dispose in a puncture resistant container. Refer to local regulatory guidelines.
<b>Symbols used:</b>	
	Expiration Date
	Lot number
	Reference number
	electronic Instruction For Use, see under <a href="http://www.gemma-medical.ch">www.gemma-medical.ch</a>
	CE-marked
	Do not reuse
	Do not use in case of broken package
	Store dry
	Keep away from sunlight
	Sterilized using ethylene oxide